





NEW PRODUCTS



Micro Vascular Profunda Clamp

DeBakey Jaws



NEW

SCANLAN®
Micro Vascular Clamp
Page 2.6

| Profunda Clamp | |
|----------------|------------------|
| 3003-617 | 5-1/4" (13.5 cm) |

Micro Vascular Lambert-Kay Clamp

Angled Shanks

Cooley Jaws



NEW

SCANLAN®
Micro Vascular Clamp
Page 2.7

| Lambert-Kay Clamp | |
|-------------------|----------------|
| 3003-59 | 4-1/2" (11 cm) |



LEGACY Forceps with SCAN-II® Handle

The SCAN-II® Handle is a revolutionary design available in the LEGACY series that features a flat surface for control with round knurled sides for a secured grip and the ability to roll the instrument during procedures.



SCANLAN® LEGACY Forceps with SCAN-II® Handle

Platform Ring Tip
 Page 3.5



| 1.0 mm Straight | |
|-----------------|------------------|
| 3003-942 | 7-1/8" (18.1 cm) |
| 3003-944 | 8-1/4" (21 cm) |
| 3003-946 | 9" (23 cm) |
| 3003-948 | 10" (25.4 cm) |

**Platform
 DIAMOND DUST™ Tip**
 Page 3.5



| 0.5 mm Straight | |
|-----------------|------------------|
| 3003-962 | 7-1/8" (18.1 cm) |
| 3003-964 | 8-1/4" (21 cm) |
| 3003-966 | 9" (23 cm) |
| 3003-968 | 10" (25.4 cm) |

| 1.0 mm Straight | |
|-----------------|------------------|
| 3003-972 | 7-1/8" (18.1 cm) |
| 3003-974 | 8-1/4" (21 cm) |
| 3003-976 | 9" (23 cm) |
| 3003-978 | 10" (25.4 cm) |

DeBakey Tip
 Page 3.6



| 1.0 mm Straight | |
|-----------------|------------------|
| 3003-952 | 7-1/8" (18.1 cm) |
| 3003-954 | 8-1/4" (21 cm) |
| 3003-956 | 9" (23 cm) |
| 3003-958 | 10" (25.4 cm) |

1x2 Micro Teeth Tip
 Page 3.7



| 0.5 mm Straight | |
|-----------------|------------------|
| 3003-980 | 7-1/8" (18.1 cm) |
| 3003-984 | 8-1/4" (21 cm) |
| 3003-990 | 9" (23 cm) |
| 3003-992 | 10" (25.4 cm) |

| 1.0 mm Straight | |
|-----------------|----------------|
| 3003-988 | 8-1/4" (21 cm) |
| 3003-994 | 9" (23 cm) |
| 3003-996 | 10" (25.4 cm) |

| 0.5 mm Curved | |
|-----------------|----------------|
| 3003-985 | 8-1/4" (21 cm) |

| 0.5 mm Angled | |
|-----------------|----------------|
| 3003-986 | 8-1/4" (21 cm) |

| 1.0 mm Angled | |
|-----------------|----------------|
| 3003-989 | 8-1/4" (21 cm) |

LEGACY Bayonet Style Forceps with SCAN-II® Handle

The SCAN-II® Handle is a revolutionary design available in the LEGACY series that features a flat surface for control with round knurled sides for a secured grip and the ability to roll the instrument during procedures.



1x2 Micro Teeth Tip
 Page 3.8



| | |
|-----------------|----------------|
| 0.5 mm Straight | |
| 3003-820 | 8-3/4" (22 cm) |

DeBakey Forceps with 1x2 Teeth



1x2 Multi Tooth Tip
 Page 3.16



| | |
|-----------------|------------------|
| 1.2 mm Straight | |
| 3003-051 | 7-3/4" (19.5 cm) |
| 3003-071 | 9-1/2" (24 cm) |



Super-Micro Fine Forceps

SCANLAN® Super-Micro Fine Forceps are designed with a low profile tapered shank for improved visibility under the microscope. The high quality titanium makes the instrument light weight, durable and non-magnetic to needles. The 0.3 mm DIAMOND DUST™ tip is ideal for picking up 9-0 suture (100 micron needles).



SCANLAN®
Super-Micro Fine Forceps
Round and Flat Handles



Round Handle
Page 3.10



0.3 mm Straight

| | |
|-----------------|------------------|
| 3303-470 | 7-1/4" (18.5 cm) |
|-----------------|------------------|

0.3 mm Angled

| | |
|-----------------|--------------|
| 3303-452 | 5" (12.7 cm) |
|-----------------|--------------|

0.3 mm Curved

| | |
|-----------------|------------------|
| 3303-477 | 7-1/4" (18.5 cm) |
|-----------------|------------------|



Flat Handle
Page 3.15



0.3 mm Straight

| | |
|-----------------|--------------|
| 3303-410 | 5" (12.7 cm) |
|-----------------|--------------|

| | |
|-----------------|------------------|
| 3303-430 | 7-1/4" (18.5 cm) |
|-----------------|------------------|

0.3 mm Angled

| | |
|-----------------|------------------|
| 3303-434 | 7-1/4" (18.5 cm) |
|-----------------|------------------|

LEGACY Needle Holders with SCAN-II® Handle

The SCAN-II® Handle is a revolutionary design available in the LEGACY series that features a flat surface for control with round knurled sides for a secured grip and the ability to roll the instrument during procedures.



SCANLAN® LEGACY Needle Holders with SCAN-II® Handle

Page 4.5



| 7-0 Suture Size | |
|-----------------|---------------|
| 3003-372 | 7" (18 cm) |
| 3003-374 | 8" (20 cm) |
| 3003-376 | 9" (23 cm) |
| 3003-378 | 10" (25.4 cm) |

Page 4.6



| 5-0 Suture Size | |
|-----------------|---------------|
| 3003-382 | 7" (18 cm) |
| 3003-384 | 8" (20 cm) |
| 3003-386 | 9" (23 cm) |
| 3003-388 | 10" (25.4 cm) |

Page 4.7



| 3-0 Suture Size | |
|-----------------|---------------|
| 3003-332 | 7" (18 cm) |
| 3003-334 | 8" (20 cm) |
| 3003-336 | 9" (23 cm) |
| 3003-338 | 10" (25.4 cm) |

NEW

Super-Micro Fine Needle Holders

SCANLAN® Super-Micro Fine Needle Holders are designed with a low profile tapered shank for improved visibility under the microscope. The high quality titanium makes the instrument light weight, durable and non-magnetic to needles. The 0.3 mm DIAMOND DUST™ tip is ideal for picking up 9-0 suture (100 micron needles).



SCANLAN®
Super-Micro Fine Needle Holders
 Round and Flat Handles

Round Handle
 Page 4.10

Flat Handle
 Page 4.14

Curved

Straight

Straight



| | |
|-----------------|--------------|
| 9-0 Suture | |
| 3303-646 | 5" (12.7 cm) |

| | |
|-----------------|------------------|
| 9-0 Suture | |
| 3303-664 | 7-1/4" (18.5 cm) |

| | |
|-----------------|------------------|
| 9-0 Suture | |
| 3303-622 | 7-1/4" (18.5 cm) |

8 inch (20 cm)

Premier™ Spring Style Scissors

The Scanlan® Premier™ spring style scissors use a revolutionary, Scanlan proprietary technology that enhances the hardness and corrosion resistant properties of stainless steel. As a result, blades retain their sharpness longer and are very resistant to environmental damage.



Round Handle
Page 5.2

| 25° Nano Blades | 45° Nano Blades | 60° Nano Blades | 90° Nano Blades |
|-----------------|-----------------|-----------------|-----------------|
| 7007-590-8 | 7007-592-8 | 7007-594-8 | 7007-596-8 |

Flat Handle
Page 5.12

| 25° Nano Blades | 45° Nano Blades | 60° Nano Blades | 90° Nano Blades |
|-----------------|-----------------|-----------------|-----------------|
| 7007-540-8 | 7007-542-8 | 7007-544-8 | 7007-546-8 |

Premier™ Super-Micro Fine Spring Style Scissors

The Scanlan® Premier™ Super-Micro Fine spring style scissors are ideal for microsurgery and pediatric Surgeries



Flat Handle
Page 5.10

| Curved Nano Blades | |
|--------------------|------------------|
| 7007-712 | 5" (12.7 cm) |
| 7007-729 | 7-1/4" (18.5 cm) |



M.D.™ Retractor System

The SCANLAN® M.D.® Retractor System featuring the new Scan-Vu® single-use radiolucent arms... Perfect for procedures performed in the hybrid operating rooms.



The SCANLAN® M.D.® Retractor System with the reusable M.D.® titanium retrator arms, are designed for light weight, reliable retraction for open surgical procedures.



SCANLAN®
M.D.® Retractor System
Page 6.20

NEW

NEW PRODUCTS

Gening J SPH Snake™ Dissector

The SCANLAN® Gening J SPH Snake™ Dissector has triangular shaped tips for precise control during dissection and the size and shape of the tip are optimized for lymph node grasping.



SCANLAN®
Gening J SPH Snake™ Dissector
Page 7.15

| With Ratchet | | Without Ratchet | |
|--------------|------------------|-----------------|------------------|
| 9009-284 | 13 1/2" (34 cm) | 9009-285 | 13 1/2" (34 cm) |
| 9009-286 | 16 3/4 (42.4 cm) | 9009-287 | 16 3/4 (42.4 cm) |




Thoroscopic Grasping Dissector

NEW

The ratchets on this instrument are tuned to engage after the jaws meet making it function like a non-ratcheted dissector. The oval rings provide a secure grasp of lymph nodes during dissection.

The curved jaw improves visualization during VATS procedures.



SCANLAN®
Thoroscopic Grasping Dissector
Page 7.12

With Ratchet

| | |
|----------|-------------|
| 9009-290 | 13" (33 cm) |
|----------|-------------|

Periosteal Elevator

NEW

Spatula like instrument for removing tissue from rib prior to resection.



SCANLAN®
Periosteal Elevator
Page 7.46


| | |
|-------------|---------------|
| 9009-823-14 | 14" (35.6 cm) |
|-------------|---------------|

| | |
|-------------|---------------|
| 9009-823-17 | 17" (43.2 cm) |
|-------------|---------------|

Foerster Grasping Clamp

NEW

The 25 x 35 oval ring jaws provide an atraumatic, secure grasp on delicate, friable tissue such as emphysematous tissue. The classic length works well for a variety of patient sizes and for most multi-port VATS procedures.



SCANLAN®
Foerster Lung Grasping Clamp
Page 7.10

Without Ratchet

| | |
|------------|-------------------|
| 9909-920-1 | 13-1/4" (33.5 cm) |
|------------|-------------------|

NEW

Scanturian® Minimally Invasive Surgical Instruments

The SCANLAN® Scanturian® MIS Instruments are individually handcrafted assuring superior quality, durability and performance you expect in your Scanlan instruments.

— FORCEPS —

DeBakey Jaws | Resano Jaws
Page 8.1



— NEEDLE HOLDERS —

Carbon Inlay Jaws
Page 8.2



— NEEDLE HOLDERS —

DIAMOND DUST™ Jaws
Page 8.3



— SCISSORS —

Straight | Curved | Angled
Page 8.4

