

TITANIUM SPRING STYLE NEEDLE HOLDERS

LEGACY Needle Holders


LEGACY Round Handle 4.1
 LEGACY Flat Handle 4.5

Titanium Round Handle

Super-Micro Fine..... 4.8
 Micro Jaws 4.9
 Heavy Jaws..... 4.13

Titanium Flat Handle

Super-Micro Fine..... 4.14
 Micro Jaws 4.15
 Heavy Jaws..... 4.18



TITANIUM RING HANDLE

Micro Jaws 4.19
 Heavy Jaws 4.25

TITANIUM BAYONET STYLE

Micro Jaws 4.26

SPRING STYLE NEEDLE HOLDERS

Round Handle

Micro Jaws 4.27
 Heavy Jaws 4.34

Flat Handle

Micro Jaws 4.35

RING HANDLE NEEDLE HOLDERS

Micro Jaws 4.39
 Fine Jaws 4.41
 Heavy Jaws 4.45

BAYONET STYLE

Micro Jaws 4.47



Needle Holder Modifications

Modification	Description
905	Cut down jaws/tips
905E	Cut down jaws/tips 0.4 mm
905TC	Cut down carbon-inlay jaws/tip
915	Gold plating
923	Satin finish instrument
924	Tension modification
928	Modify to left-handed instrument (ring handle)
928A	Modify to left-handed instrument (spring style)
929	Special etching
971	Lengthen instrument
973	Shorten instrument




3003-204

3003-204 with modifications **971**



3003-180

3003-180 with modifications **905**



6006-32

6006-32 with modifications **928**

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.




Features of *LEGACY* Titanium Needle Holders

Titanium NEEDLE HOLDERS

Unique 10 mm Diameter and Extended Knurling

- Provides the titanium instruments with the optimized weight and balance.
- The extended knurling provides more surface area for the surgeon to grasp the instrument.




Perfect-Pyramid™

Provides a secure grasp.

Non-Magnetic

The properties of titanium prevent the needles from clinging to the instrument.



SCAN-II® Handle

Revolutionary design that features a flat surface for control with round knurled sides for a secured grip and the ability to roll the instrument during procedures.

Diamond-Dust™

Securely and delicately grips the needles to prevent slippage and any damage to the needle.



Round Handle

The original 10 mm round design provides the surgeon with a secure grip and the ability to roll the instrument during procedures.



3003-342

Instrument is shown approximate size




SCANLAN®

LEGACY Needle Holders

Round Handle | 10 mm Diameter, Extended Knurling

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-342	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	7" (18 cm)
	3003-344			8" (20 cm)
	3003-346			9" (23 cm)
	3003-348			10" (25.4 cm)

3003-381




SCANLAN®

LEGACY Puskas™ Jacobson Micro Needle Holders

Round Handle | 8 mm Diameter, Extended Knurling

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-381	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	8 1/4" (21 cm)
	3003-391			9 1/2" (24 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Spring Style | LEGACY | Round Handle | 5-0 Suture

3003-352
Instrument is shown approximate size



SCANLAN®
LEGACY Needle Holders
Round Handle | 10 mm Diameter, Extended Knurling

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-352			7" (18 cm)
	3003-354	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	5-0 and Smaller	8" (20 cm)
	3003-356			9" (23 cm)
	3003-358			10" (25.4 cm)

3003-390



SCANLAN®
LEGACY Puskas™ Jacobson Micro Needle Holders
Round Handle | 8 mm Diameter, Extended Knurling

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-390	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	5-0 and Smaller	9" (23 cm)

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
For a complete listing of models visit... www.scanlaninternational.com





SCANLAN®

LEGACY Jacobson Micro Needle Holders

Round Handle | 10 mm Diameter, Extended Knurling

Jaws

(enlarged to show detail)



Cat. No.	Description	Recommended Suture Size	Total Length
3003-362			7" (18 cm)
3003-364	Straight DIAMOND DUST™ Jaws	3-0 and Smaller	8" (20 cm)
3003-366	Streamline Box Lock With Ratchet		9" (23 cm)
3003-368			10" (25.4 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Spring Style | LEGACY | Flat Handle | 7-0 Suture

3003-374



Titanium NEEDLE HOLDERS


SCANLAN®

LEGACY Needle Holders

SCAN-II® Handle | 10 mm, Extended Knurling

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	NEW 3003-372			7" (18 cm)
	NEW 3003-374	Straight DIAMOND DUST™ Jaws	7-0 and Smaller	8" (20 cm)
	NEW 3003-376	Regular Box Lock With Ratchet		9" (23 cm)
	NEW 3003-378			10" (25.4 cm)


Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.

For a complete listing of models visit... www.scanlaninternational.com





SCANLAN®
LEGACY Needle Holders
 SCAN-II® Handle | 10 mm, Extended Knurling

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	NEW 3003-382			7" (18 cm)
	NEW 3003-384	Straight DIAMOND DUST™ Jaws	5-0 and Smaller	8" (20 cm)
	NEW 3003-386	Streamline Box Lock With Ratchet		9" (23 cm)
	NEW 3003-388			10" (25.4 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Spring Style | LEGACY | Flat Handle | 3-0 Suture

3003-334

Titanium NEEDLE HOLDERS



SCANLAN®
LEGACY Needle Holders
 SCAN-II® Handle | 10 mm, Extended Knurling

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	NEW 3003-332	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	3-0 and Smaller	7" (18 cm)
	NEW 3003-334			8" (20 cm)
	NEW 3003-336			9" (23 cm)
	NEW 3003-338			10" (25.4 cm)

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
 For a complete listing of models visit... www.scanlaninternational.com



3303-646

Instrument is shown approximate size



SCANLAN®
Super-Micro Fine Needle Holders
 Round Handle

Jaws
 (enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	NEW 3303-646	Curved 0.3 mm DIAMOND DUST™ Jaws	9-0 and Smaller	5" (12.7 cm)
	3303-667	Regular Box Lock Without Ratchet		7 1/4" (18.5 cm)

3303-664



SCANLAN®
Super-Micro Fine Needle Holders
 Round Handle

Jaws
 (enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	NEW 3303-664	Straight 0.3 mm DIAMOND DUST™ Jaws	9-0 and Smaller	7 1/4" (18.5 cm)
		Regular Box Lock Without Ratchet		

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Spring Style | Round Handle | 8-0 Suture



Titanium NEEDLE HOLDERS


SCANLAN®

Jacobson Micro Needle Holders

Round Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-177-3	0.3 mm Tapered DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	8-0 and Smaller	7 1/4" (18.5 cm)
	3003-177-4	0.4 mm Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.

For a complete listing of models visit... www.scanlaninternational.com





SCANLAN®
Jacobson Micro Needle Holders

Round Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-171	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	7 1/4" (18.5 cm)
	3003-173	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	3003-175	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	3003-177	Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	3003-181	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		8" (20 cm)
	3003-183	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	3003-185	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	3003-187	Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	3003-191	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		9" (23 cm)
	3003-193	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	3003-195	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		

Spring Style | Round Handle | 5-0 Suture

3003-93

Instrument is shown approximate size



Titanium NEEDLE HOLDERS

SCANLAN®
Castroviejo Micro Needle Holders
 Round Handle

Jaws
 (enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-90	Straight DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet	5-0 and Smaller	4 3/4" (12 cm)
	3003-90A	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		
	3003-91	Curved DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		
	3003-91A	Curved DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		
	3003-93	Straight DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		6" (15 cm)
	3003-93A	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		
	3003-94	Curved DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		
	3003-94A	Curved DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		
	3003-94-4	Curved 0.4 mm DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		
	3003-94-4	Curved 0.4 mm DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		

Please Note...
 Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible.
 Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
 For a complete listing of models visit... www.scanlaninternational.com





SCANLAN®

Jacobson Needle Holders

Round Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-170	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	5-0 and Smaller	7 1/4" (18.5 cm)
	3003-172	Straight DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		
	3003-174	Curved DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		
	3003-180	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		
	3003-182	Straight DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		8" (20 cm)
	3003-184	Curved DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		
	3003-190	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		9" (23 cm)
	3003-200	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		10" (25.4 cm)
	3003-202	Straight DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		
	3003-204	Curved DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		
	3003-200-12	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	12" (30.5 cm)	

Spring Style | Round Handle | 2-0 Suture


3003-530

Titanium NEEDLE HOLDERS



SCANLAN®
Jacobson Needle Holders
 For Large Needles | Round Handle

Jaws
 (enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-510			7 1/2" (19 cm)
	3003-530	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	2-0 and Smaller	9" (23 cm)
	3003-540			10" (25.4 cm)
	NEW 3003-540-12			12" (30.5 cm)

Please Note...
 Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
 For a complete listing of models visit... www.scanlaninternational.com



3303-603

Instrument is shown approximate size



3303-622




SCANLAN®

Super-Micro Fine Needle Holders

Round Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3303-603	Straight 0.3 mm DIAMOND DUST™ Jaws	9-0 and Smaller	5" (12.7 cm)
	NEW 3303-622	Regular Box Lock Without Ratchet		7 1/4" (18.5 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

3003-225

Titanium NEEDLE HOLDERS






SCANLAN®

Jacobson Micro Needle Holders

Flat Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-207	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	5 3/4" (14.5 cm)
	3003-225	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		7 1/4" (18.5 cm)
	3003-227	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		8" (20 cm)
	3003-229	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	3003-261	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	3003-263	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	3003-265	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.

For a complete listing of models visit... www.scanlaninternational.com










SCANLAN®

Jacobson Micro Needle Holders

Flat Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-235	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	8 1/4" (21 cm)
	3003-237	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	3003-239	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	3003-234-1	Straight Fine DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		8 1/2" (21.5 cm)
	3003-243	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		9 1/2" (24 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

3003-224



Titanium NEEDLE HOLDERS






SCANLAN®

Jacobson Micro Needle Holders

Flat Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-224	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	5-0 and Smaller	7 1/4" (18.5 cm)
	3003-226	Straight DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		
	3003-234	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		
	3003-236	Straight DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		8 1/2" (21.5 cm)
	3003-238	Curved DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.

For a complete listing of models visit... www.scanlaninternational.com





SCANLAN®

Rendina Needle Holders

For Large Needles | Flat Handle Anniversary Knurl

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-550			7 1/2" (19 cm)
	3003-570	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	2-0 and Smaller	9" (23 cm)
	3003-580			10" (25.4 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

3003-450

Instrument is shown approximate size



Titanium NEEDLE HOLDERS

SCANLAN®
Richards Micro Needle Holders
 Ring Handle

Jaws
 (enlarged to show detail)


	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-450			6" (15 cm)
	3003-452	Straight Extra Fine DIAMOND DUST™ Jaws	6-0 and Smaller	6 3/4" (17 cm)
	3003-454			7 1/2" (19 cm)

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
 For a complete listing of models visit... www.scanlaninternational.com



3003-430


Instrument is shown approximate size



SCANLAN®
Ryder Needle Holders
 Ring Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-430			6" (15 cm)
	3003-432	Straight DIAMOND DUST™ Jaws	5-0 and Smaller	6 3/4" (17 cm)
	3003-434			7 1/2" (19 cm)

Please Note...


Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

3003-412



Titanium NEEDLE HOLDERS

SCANLAN®
DeBakey Micro Needle Holders
Ring Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-412	Straight DIAMOND DUST™ Jaws	5-0 and Smaller	7 1/8" (18.1 cm)
	3003-414			8" (20 cm)
	3003-416			8 5/8" (21.9 cm)
	3003-418			10" (25.4 cm)

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
For a complete listing of models visit... www.scanlaninternational.com





SCANLAN®

Cardiovascular Needle Holders

Ring Handle

Jaws

(enlarged to show detail)


	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-436	Straight Ryder Style DIAMOND DUST™ Jaws	5-0 and Smaller	7 1/8" (18.1 cm)
	3003-438			8" (20 cm)
	3003-440			8 5/8" (21.9 cm)
	3003-442			10" (25.4 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.


3003-422

Titanium NEEDLE HOLDERS



SCANLAN®
Microvascular Needle Holders
 Ring Handle

Jaws
 (enlarged to show detail)




	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-422	Straight DIAMOND DUST™ Jaws	5-0 and Smaller	7 1/8" (18.1 cm)
	3003-424			8" (20 cm)
	3003-426			8 5/8" (21.9 cm)
	3003-428			10" (25.4 cm)

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
 For a complete listing of models visit... www.scanlaninternational.com





SCANLAN®
Micro Needle Holders
Ring Handle


Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-462	Straight DIAMOND DUST™ Jaws	5-0 and Smaller	7 1/4" (18.5 cm)
	3003-466			8 3/4" (22 cm)
	3003-464	Curved DIAMOND DUST™ Jaws		7 1/4" (18.5 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.


3003-402

Instrument is shown approximate size



Titanium NEEDLE HOLDERS

SCANLAN®
Mayo-Hegar Needle Holders
 Ring Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-402	Straight DIAMOND DUST™ Jaws	2-0 and Smaller	6 1/4" (16 cm)
	3003-404			7 1/8" (18.1 cm)
	3003-406			8" (20 cm)
	3003-408			8 5/8" (21.9 cm)

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
 For a complete listing of models visit... www.scanlaninternational.com




Bayonet Style | Round Handle | 5-0 Suture

Titanium NEEDLE HOLDERS

3789-11T



SCANLAN®
Micro Needle Holders
Round Handle





Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	3789-11T	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	5-0 and Smaller	9 1/2" (24 cm)

Bayonet Style | Flat Handle | 5-0 Suture

3003-866



SCANLAN®
Micro Needle Holders
Flat Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	3003-866	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	6-0 and Smaller	8 3/4" (22 cm)
	3003-868	Straight DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		
	3003-870	Curved DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		
	3003-872	Curved DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		

Spring Style | Round Handle | 8-0 Suture

6006-56-3
Instrument is shown approximate size



Titanium NEEDLE HOLDERS





SCANLAN®
Castroviejo Micro Needle Holders
Round Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-56-3	Curved 0.3 mm Tapered DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	8-0 and Smaller	5 3/4" (14.5 cm)

6006-44
Instrument is shown approximate size




SCANLAN®
Castroviejo Micro Needle Holders
Round Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-44	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	4 3/4" (12 cm)
	6006-45	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	6006-46	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	6006-47	Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		

Please Note...
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible.
Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
For a complete listing of models visit... www.scanlaninternational.com









SCANLAN®
Jacobson Micro Needle Holders

Round Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-121	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	8 1/4" (21 cm)
	6006-123	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	6006-125	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	6006-127	Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Spring Style | Round Handle | 7-0 Suture




6006-128




SCANLAN®

Jacobson Micro Needle Holders

Short Knurling | Round Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-128	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	9" (23 cm)
	6006-130	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	6006-132	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		




6006-129



SCANLAN®

Jacobson Micro Needle Holders

Round Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-129	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	9" (23 cm)
	6006-131	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	6006-133	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
For a complete listing of models visit... www.scanlaninternational.com



6006-48

Instrument is shown approximate size





SCANLAN®

Glasser Micro Needle Holders

Inside Spring | Round Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-48	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	5-0 and Smaller	4 7/8" (12.4 cm)
	6006-139	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		7" (18 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.



Spring Style | Round Handle | 5-0 Suture

6006-50

Instrument is shown approximate size



Stainless Steel NEEDLE HOLDERS



SCANLAN®

Castroviejo Micro Needle Holders

Round Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-50	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	5-0 and Smaller	5 3/4" (14.5 cm)
	6006-52	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	6006-54	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	6006-56	Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.

For a complete listing of models visit... www.scanlaninternational.com



6006-59

Instrument is shown approximate size







SCANLAN®

Jacobson Micro Needle Holders

Round Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-59	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	5-0 and Smaller	6 1/4" (16 cm)
	6006-61	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	6006-63	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	6006-65	Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		

Please Note...





Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Spring Style | Round Handle | 5-0 Suture

6006-82






SCANLAN®
Jacobson Micro Needle Holders
 Round Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-82	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	5-0 and Smaller	7 1/4" (18.5 cm)
	6006-84	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	6006-86	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	6006-88	Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		

6006-120




SCANLAN®
Jacobson Micro Needle Holders
 Round Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-120	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	5-0 and Smaller	8 1/4" (21 cm)
	6006-122	Straight DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet		
	6006-124	Curved DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
 For a complete listing of models visit... www.scanlaninternational.com







SCANLAN®

Jacobson Heavy Jaw Needle Holders

For Large Needles | Round Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-310	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	2-0 and Smaller	7 1/2" (19 cm)
	6006-330			9" (23 cm)
	6006-340			10" (25.4 cm)
	6006-334	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		9" (23 cm)


Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Spring Style | Flat Handle | 7-0 Suture

6006-70

Instrument is shown approximate size



SCANLAN®

Jacobson Micro Needle Holders

Flat Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-70	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	6 3/4" (17 cm)
	6006-72	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		

Stainless Steel NEEDLE HOLDERS

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.

For a complete listing of models visit... www.scanlaninternational.com



6006-58

Instrument is shown approximate size





SCANLAN®

Castroviejo Micro Needle Holders

Flat Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-58	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	5-0 and Smaller	5 1/2" (14 cm)
	6006-62	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		

6006-90






SCANLAN®

Jacobson Micro Needle Holders

Flat Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-90	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	5-0 and Smaller	7 1/4" (18.5 cm)
	6006-92	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	6006-94	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Spring Style | Flat Handle | 5-0 Suture

6006-108





SCANLAN®

Jacobson Micro Needle Holders

Flat Handle

Jaws


(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-108	Straight DIAMOND DUST™ Jaws	5-0 and Smaller	8 1/4" (21 cm)
	6006-158	Regular Box Lock With Ratchet		9 1/2" (24 cm)
	6006-110	Straight DIAMOND DUST™ Jaws		8 1/4" (21 cm)
	6006-159	Regular Box Lock Without Ratchet		9 1/2" (24 cm)
	6006-112	Curved DIAMOND DUST™ Jaws		8 1/4" (21 cm)
	6006-160	Regular Box Lock With Ratchet		9 1/2" (24 cm)

Stainless Steel NEEDLE HOLDERS

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
For a complete listing of models visit... www.scanlaninternational.com







SCANLAN®

Jacobson Micro Needle Holders

Flat Knurled Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-170		5-0 and Smaller	7 1/4" (18.5 cm)
	6006-100	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		8" (20 cm)
	6006-174			9" (23 cm)
	6006-102	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		8" (20 cm)
	6006-104	Curved DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		8" (20 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

4157-77

Instrument is shown approximate size



SCANLAN®
Kaye Needle Holders / Scissors
Ring Handle

Jaws
(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	4157-77	Straight DIAMOND DUST™ Jaws	5-0 and Smaller	4 3/4" (12 cm)

4155-07

Instrument is shown approximate size



SCANLAN®
Webster Micro Needle Holders
Ring Handle

Jaws
(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	4155-06	Straight Delicate DIAMOND DUST™ Jaws	5-0 and Smaller	5 1/8" (13 cm)
	4155-07	Straight Delicate DIAMOND DUST™ Jaws Offset Bows		

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
For a complete listing of models visit... www.scanlaninternational.com





SCANLAN®

Jacobson Micro Needle Holders

Ring Handle

Jaws

(enlarged to show detail)


	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-148	Straight DIAMOND DUST™ Jaws With Ratchet	5-0 and Smaller	7 1/4" (18.5 cm)
	6006-152			9" (23 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

6006-04

Stainless Steel NEEDLE HOLDERS



SCANLAN®
Dennis Delicate Needle Holders
 Ring Handle

Jaws <small>(enlarged to show detail)</small>	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-04	Straight DIAMOND DUST™ Jaws	5-0 and Smaller	7" (18 cm)



Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
 For a complete listing of models visit... www.scanlaninternational.com





SCANLAN®
Micro Needle Holders
 Ryder-Style | Ring Handle


Jaws
 (enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-02	Straight 1.5 mm DIAMOND DUST™ Jaws	5-0 and Smaller	6 1/4" (16 cm)
	6006-29	Straight 1.2 mm DIAMOND DUST™ Jaws		6 1/4" (16 cm)
	6006-33			7 1/4" (18.5 cm)
	6006-37			8 3/4" (22 cm)
	6006-41			10" (25.4 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.


6006-32



Stainless Steel NEEDLE HOLDERS

SCANLAN®
Micro Needle Holders
 Ring Handle

Jaws
 (enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-20	Straight DIAMOND DUST™ Jaws	5-0 and Smaller	5" (12.7 cm)
	6006-24			5 3/4" (14.5 cm)
	6006-28			6 1/4" (16 cm)
	6006-32			7 1/4" (18.5 cm)
	6006-36			8 3/4" (22 cm)
	6006-40			10" (25.4 cm)

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
 For a complete listing of models visit... www.scanlaninternational.com



6006-08

Instrument is shown approximate size



SCANLAN®
Fell Needle Holders
 Ring Handle

Jaws
 (enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-08	2 mm Straight Carbon-Inlay Jaws Gold Handle	3-0 and Smaller	6" (15 cm)
	6006-10			7" (18 cm)
	6006-12			8" (20 cm)
	6006-14			8 1/2" (21.5 cm)
	6006-17			10 1/2" (27 cm)
	6006-19			11 3/4" (29.8 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

6006-630

Stainless Steel NEEDLE HOLDERS




SCANLAN®

Heavy Jaw Needle Holders

Ring Handle

Jaws

(enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-610			7 1/4" (18.5 cm)
	6006-620	Straight DIAMOND DUST™ Jaws	2-0 and Smaller	8 1/8" (20.1 cm)
	6006-630			9 1/8" (23.5 cm)

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.

For a complete listing of models visit... www.scanlaninternational.com





Ring Handle | 2-0 Suture

Stainless Steel NEEDLE HOLDERS

6006-16



SCANLAN®
Stratte Valve Needle Holders
Ring Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	6006-15	Curved Carbon-Inlay Jaws Gold Handle	2-0 and Smaller	7 1/2" (19 cm)
	6006-16			10" (25.4 cm)

Ring Handle | #5 Wire

6006-05



SCANLAN®
Sternal Needle Holders
Ring Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Wire Size	Total Length
	6006-05	Straight Carbon-Inlay Jaws Gold Handle	#5 and Smaller	7" (18 cm)

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.

For a complete listing of models visit... www.scanlaninternational.com

Bayonet Style | Round Handle | 8-0 Suture

3789-14-3



SCANLAN®
Micro Needle Holders
 Round Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	3789-14-3	Curved Up 0.3 mm DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	8-0 and Smaller	9" (23 cm)





SCANLAN®

Bayonet Style | Round Handle | 7-0 Suture

3789-11



SCANLAN®
Micro Needle Holders
 Round Handle

Jaws (enlarged to show detail)	Cat. No.	Description	Recommended Suture Size	Total Length
	3789-11	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	9" (23 cm)
	3789-13	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	3789-12	Curved Up DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	3789-14	Curved Up DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		





Stainless Steel NEEDLE HOLDERS





SCANLAN®
Micro Needle Holders
 Flat Handle

Jaws
 (enlarged to show detail)

	Cat. No.	Description	Recommended Suture Size	Total Length
	3788-12	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	9" (23 cm)
	3788-14	Straight DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		
	3788-13	Curved Up DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		
	3788-15	Curved Up DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet		

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.

For a complete listing of models visit... www.scanlaninternational.com