



Scanturian® Forceps ..... 8.1

Scanturian® Needle Holders ..... 8.2

Scanturian® Scissors ..... 8.4

Scanturian® V.I.P.™ Sterilization Trays ..... 8.6

V.I.P.™ Cleaning Tray ..... 8.8



# Scanturian® Minimally Invasive Surgical Instruments

MIS



### ScanShield® Shaft Feature

The ScanShield® shaft feature surrounds the hinge mechanisms to avoid suture entrapment and intraoperative trauma to valve tissue and chordae tendineae (patent pending)



### Tapered Shaft

The tapered shaft design allows for more robust moving mechanisms and efficient cleaning due to the increased shaft diameter at the distal tip



### Scanturian® MIS Scissors with SUPER CUT™ Blades

A unique blade design - combining a razor edge with a scissors edge, these scissors produce a clean, forward cut, eliminating pinching or slipping of tissue

### Instrument Identification

Color-coded and etching for easy intraoperative instrument identification

#### FORCEPS



#### NEEDLE HOLDERS



#### NEEDLE HOLDERS



#### SCISSORS



### Titanium Handle

Well balanced, durable and handcrafted MIS Instrument design

### Handle Knurling

Proven handle knurling design provides comfortable, secure grip and ergonomic handle design

### Proximal End Cleaning Port

Luer lock flush port system allows for reliable automated and intraoperative instrument flushing

### Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.



# Scanturian® Forceps

MIS

9019-420-17



9019-420-25



Cat. No.	Description	Distal Tip	Working Length	Total Length	Handle Color
9019-420-17	SCANLAN® Scanturian® MIS DeBakey Forceps 1.5 mm Straight Jaws Without Ratchets		6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-420-25			9 7/8" (25 cm)	15 1/3" (39 cm)	
9019-430-17	SCANLAN® Scanturian® MIS DeBakey Forceps 2.5 mm Straight Jaws Without Ratchets		6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-430-25			9 7/8" (25 cm)	15 1/3" (39 cm)	
9019-440-17	SCANLAN® Scanturian® MIS Resano Forceps 2.5 mm Straight Jaws Without Ratchets		6 3/4" (17 cm)	12- 1/4" (31 cm)	
9019-440-25			9 7/8" (25 cm)	15 1/3" (39 cm)	

Please Note...  
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible.  
Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.  
For a complete listing of models visit... [www.scanlaninternational.com](http://www.scanlaninternational.com)



# Scanturian® Needle Holders

MIS

9019-610-17



9019-610-25



**MIS** SCANLAN®  
Scanturian® Needle Holders

Cat. No.	Description	Distal Tip	Suture Size	Working Length	Total Length	Handle Color	
9019-610-17	SCANLAN® Scanturian® MIS Needle Holder Curved Carbon-Inlay Jaws with Ratchet		2-0 and Smaller	6 3/4" (17 cm)	12 1/4" (31 cm)		
9019-610-25				9 7/8" (25 cm)	15 1/3" (39 cm)		
9019-614-17	SCANLAN® Scanturian® MIS Needle Holder Straight Carbon-Inlay Jaws with Ratchet			6 3/4" (17 cm)	12 1/4" (31 cm)		
9019-614-25				9 7/8" (25 cm)	15 1/3" (39 cm)		
9019-630-17	SCANLAN® Scanturian® MIS Needle Holder Straight Carbon-Inlay Jaws with Ratchet			4-0 and Smaller	6 3/4" (17 cm)		12 1/4" (31 cm)
9019-630-25				9 7/8" (25 cm)	15 1/3" (39 cm)		

**Please Note...**

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

# Scanturian® Needle Holders

9019-640-17




9019-640-25



MIS

**MIS** SCANLAN®  
Scanturian® Needle Holders

Cat. No.	Description	Distal Tip	Suture Size	Working Length	Total Length	Handle Color
9019-640-17	SCANLAN® Scanturian® MIS Needle Holder Straight DIAMOND DUST™ Jaws with Ratchet		5-0 and Smaller	6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-640-25				9 7/8" (25 cm)	15 1/3" (39 cm)	

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.  
For a complete listing of models visit... [www.scanlaninternational.com](http://www.scanlaninternational.com)



MIS

9019-722-17SC



9019-722-25SC



**MIS** SCANLAN®  
Scanturian® Scissors

Cat. No.	Description	Distal Tip	Working Length	Total Length	Handle Color
9019-710-17SC	SCANLAN® Scanturian® MIS Super Cut™ Scissors Straight Fine Blades, 1 cm, Blunt Tips		6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-710-25SC			9 7/8" (25 cm)	15 1/3" (39 cm)	
9019-715-17SC	SCANLAN® Scanturian® MIS Super Cut™ Scissors Curved Fine Blades, 1 cm, Blunt Tips		6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-715-25SC			9 7/8" (25 cm)	15 1/3" (39 cm)	
9019-720-17SC	SCANLAN® Scanturian® MIS Super Cut™ Scissors 15° Curved Blades		6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-720-25SC			9 7/8" (25 cm)	15 1/3" (39 cm)	
9019-722-17SC	SCANLAN® Scanturian® MIS Super Cut™ Scissors 30° Curved Blades		6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-722-25SC			9 7/8" (25 cm)	15 1/3" (39 cm)	
9019-724-17SC	SCANLAN® Scanturian® MIS Super Cut™ Scissors 70° Curved Blades		6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-724-25SC			9 7/8" (25 cm)	15 1/3" (39 cm)	

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

SCANLAN®  
**Scanturian® Scissors**

MIS





9019-740-17SC



9019-740-25SC



MIS SCANLAN®  
**Scanturian® Scissors**

Cat. No.	Description	Distal Tip	Working Length	Total Length	Handle Color
9019-740-17SC	SCANLAN® Scanturian® Potts Scissors 45° Angled Fine Blades		6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-740-25SC			9 7/8" (25 cm)	15 1/3" (39 cm)	
9019-748-17SC	SCANLAN® Scanturian® Potts Scissors 90° Angled Fine Blades		6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-748-25SC			9 7/8" (25 cm)	15 1/3" (39 cm)	
9019-750-17SC	SCANLAN® Scanturian® Potts Scissors 125° Angled Fine Blades		6 3/4" (17 cm)	12 1/4" (31 cm)	
9019-750-25SC			9 7/8" (25 cm)	15 1/3" (39 cm)	

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.  
 For a complete listing of models visit... [www.scanlaninternational.com](http://www.scanlaninternational.com)





SCANLAN®

## V.I.P.™ Sterilization Trays

Metal Sterilization Trays, Aluminum, Scan-Guard™ Silicone Inserts, Detachable Lids


Cat. No.	Description	Dimensions	Capacity
2160-03	SCANLAN® Scanturian® Valued Instrument Protection (V.I.P.™) Sterilization Tray	21 3/4" x 10 1/4" x 2 1/2" (55.2cm x 26cm x 6.4cm)	Up to 10 MIS Instruments

All trays are steam, Immediate Use Steam Sterilization (IUSS) and ETO sterilization compatible  
 Firm Sale - May be available for evaluation

Please Note...

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.





SCANLAN®

**V.I.P.™ Sterilization Trays**

Metal Sterilization Trays, Aluminum, Scan-Guard™ Silicone Inserts, Detachable Lids

Cat. No.	Description	Dimensions	Capacity
2161-03	SCANLAN® Scanturian® Valued Instrument Protection (V.I.P.™) Sterilization Tray with additional tray insert	21 3/4" x 10 1/4" x 4" (55.2 cm x 26 cm x 10 cm)	Up to 18 MIS Instruments

All trays are steam, Immediate Use Steam Sterilization (IUSS) and ETO sterilization compatible  
Firm Sale - May be available for evaluation

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.  
For a complete listing of models visit... [www.scanlaninternational.com](http://www.scanlaninternational.com)





SCANLAN®  
**V.I.P.™ Cleaning Tray**  
 Metal Cleaning Tray, Aluminum, Scan-Guard™ Silicone Inserts

Cat. No.	Description	Dimensions	Capacity
2162-03	SCANLAN® Valued Instrument Protection (V.I.P.™) Cleaning Tray	18 7/8" x 9 7/8" x 2 1/4" (48cm x 25cm x 5.8cm)	Up to 10 instruments

Firm Sale - May be available for evaluation

**Please Note...**

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.  
 For a complete listing of models visit... [www.scanlaninternational.com](http://www.scanlaninternational.com)